

Voltammetric determination of mexidol

Ziyatdinova G., Budnikov G.

Kazan Federal University, 420008, Kremlevskaya 18, Kazan, Russia

Abstract

A new voltammetric method for determining mexidol in pharmaceuticals is proposed. Using this technique, mexidol was determined in model solutions (with RS = of 1-6%) and in a ready-to-use preparation. The analytical range of mexidol determination using a glassy carbon composite electrode in 0.1 M H₂SO₄ extends from 4.8×10^{-3} to 1.8×10^{-2} M. The detection limit is 1.9×10^{-3} M. The proposed procedure can be used for the quality control in the drug production technology. © 2005 Springer Science+Business Media, Inc.
